# Nox Nasal Pressure Cannula with Filter

INSTRUCTIONS FOR USE

## Table of Contents

EN - English .....	3
CS - Čeština .....	5
DA - Dansk .....	7
DE - Deutsch .....	9
EL - Ελληνικά.....	11
ET - Eesti .....	14
ES - Español.....	16
FI - Suomi .....	18
FR - Français.....	20
IT - Italiano .....	22
JA - 日本語 .....	24
KO - 한국어 .....	26
LT - Lietuvių .....	28
LV - Latviski .....	30
NL - Nederlands .....	32
NO - Norsk .....	34
PL - Polski.....	36
PT - Português .....	38
SV - Svenska.....	40
ZH - 中文.....	42

### Nox Nasal Pressure Cannula with Filter – Instructions for Use

The Nox Nasal Pressure Cannula with Filter is available in one size, adult, and has a total length of 90 cm (35 in). The cannula is clean and ready for use.

Available package:

Nox Nasal Pressure Cannula with Filter, 40 units #55 20 11

#### Intended Use

The intended use of the Nox Nasal Pressure Cannula with Filter is to allow for the measuring of respiratory airflow when connected to a pressure port of a sleep recorder during a sleep study. The cannula is intended for single patient use only.

The Nox Nasal Pressure Cannula with Filter is intended for adult patients.

The Nox Nasal Pressure Cannula with Filter is intended to be used in hospitals, institutions, sleep centers, sleep clinics, or other test environments, including the patient’s home.

#### Indication for use

The Nox Nasal Pressure Cannula with Filter is an ancillary device that assists the sleep recorder in obtaining a respiratory airflow signal as one of the parameters used to aid healthcare professionals in diagnosing sleep disorders.

#### Patient Population

Patients over 18 years of age who are undergoing physiological measurements as part of a sleep test requested by a healthcare professional.

#### Intended Users

The users are healthcare professionals who have received training in the areas of hospital/clinical procedures, physiological monitoring of human subjects, or sleep disorder investigation. When conducting a sleep study in the home environment, the patient may perform a self-hook-up, with or without assistance from a family member, after being instructed by a qualified healthcare professional before being sent home with the system.

#### Clinical Benefit

The main benefit of the Nox Nasal Pressure Cannula with Filter is that it aids healthcare professionals in successfully obtaining physiological signals for a sleep test, which they can then use in the process of diagnosing sleep disorders. The device therefore provides an indirect benefit to patients through a positive impact on patient management by allowing for physiological signals to be obtained to facilitate the diagnostic process.

#### Information for Safety

The Nox Nasal Pressure Cannula with Filter is single patient use only.

#### WARNINGS:

- Do not use damaged equipment. Use of damaged equipment may result in patient harm, such as incorrect patient diagnosis.
- Re-using the cannula may result in bad quality of the signal recorded and lead to incorrect patient diagnosis.
- As with all medical equipment, carefully route cables and connections to reduce the possibility of entanglement or strangulation
- The use of the same cannula on more than one patient poses a risk of cross-contamination.
- No modification of this equipment is allowed. Any modification may lead to incorrect patient diagnosis.
- Caution: U.S. federal law restricts this device to sale by or on the order of a licensed medical practitioner.

**Note:** The Nox Nasal Pressure Cannula with Filter has a built-in hydrophobic filter and is the preferred way to measure nasal airflow and snoring as it is designed to maximize the signal quality and fits directly with the Nox Recorders.


**Note:** The Nox Nasal Pressure Cannula with Filter is validated to give high quality flow signal for use in diagnosis of sleep disorder breathing. Using a cannula that has not been validated for this purpose might compromise the signal quality and the operational reliability of the system.

#### Potential Adverse Event

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

#### Directions for Use

1. The prongs shall point downwards inside the nostrils.
2. Place the nasal prongs gently in the nostrils.
3. Pull the cannula tubing over the ears and then position it under the chin.
4. Slide the fastener snugly under the chin to hold the cannula tubing securely in place.

**Note:** Medical tape can be used to hold the cannula against the cheeks to secure the cannula in place if necessary.

The barbed end of the cannula connects to the pressure input on the Nox recorders. For further information, please refer to the applicable Nox recorders instructions for use.

#### Environmental Conditions

Temperature Operation: +5°C to +40°C (+41°F to +104°F)  
 Transportation/Storage: -25°C to 70°C (-13°F to +158°F)

### Compatible devices

Nox T3 Recorder	#56 10 10
Nox T3s Recorder	#56 10 11 / #56 11 11
Nox A1 Recorder	#56 14 10
Nox A1s Recorder	#56 14 11 / #56 11 21

(01)1569431111XXXX  
 (11)YYMMDD  
 (10)ZZZZZZ

Unique Device Identifier (UDI); Application Identifier (01) represents the device identifier (DI) (“1569431111XXXX”), the Application Identifier (11) the production date/date of manufacture (“YYMMDD”, with “YY” the last two digits of the production year, “MM” the production month and “DD” the production day), and the Application Identifier (10) the lot number of the device (“ZZZZZZ”) if applicable.

### Material Information

- Nasal prongs: Polyvinyl Chloride PVC (plastisol - liquid PVC)
- Tube: Flexible PVC (FPVC)
- Tube adapter: Flexible PVC (FPVC)
- Fastener slider: Polyethylene (PE)
- Luer-lock connector: Polymer Acrylonitrile Butadiene Styrene (ABS)
- Hydrophobic filter: Polytetrafluoroethylene (PTFE)
- Filter housing: Polypropylene (PP)
- Barbed connector: Thermoplastic Polyurethane (TPU)

### Disposal

Follow local governing ordinances and recycling instructions regarding disposal or recycling of this product. Please contact your distributor regarding take-back or recycling of the components.

The Instructions for Use are available in electronic format on Nox Medical’s website: [www.noxmedical.com/ifu](http://www.noxmedical.com/ifu). Note: A pdf reader is required to open the electronic file.

Hard copies can be requested at no additional cost by emailing [support@noxmedical.com](mailto:support@noxmedical.com). The hard copy will be sent within 7 calendar days.

### Description of Symbols and Labels

	Operating instructions / Consult instructions for use		CE marking indicating conformance to the applicable EU regulations
	Manufacturer Information		Temperature limitation
	Country of manufacture and date of manufacture		Do not re-use
	Medical device		Keep dry
	Batch code / Lot number		Fragile, handle with care
	Catalogue number / Reference number		Authorized representative for Switzerland
	Unique device identifier		Packaging unit / Number of pieces in package
	Unique Device Identifier (UDI) presented in data matrix format		Technical name
US: Rx Only	US federal law restricts this device to sale by or on the order of a licensed medical practitioner	Nox Nasal Pressure Cannula with Filter	Brand name /Model name

### Copyright Notice

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language or computer language, in any form, or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without the prior written authorization from Nox Medical.

Nox Medical ehf  
 Katrinartuni 2  
 IS-105 Reykjavík  
 Iceland  
[www.noxmedical.com](http://www.noxmedical.com)

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