

# Product Information

NOX Filter Tube Connector  
- 50 units/unités



NOX-FTC50-0200

**LOT** YYYMMDD  
Barcode

**REF**   
5 5 2 1 1 0

For distributor information go to: [www.noxmedical.com](http://www.noxmedical.com)  
Pour plus d'informations sur le distributeur : [www.noxmedical.com](http://www.noxmedical.com)

**Material Information/Information matériaux:**

Connector/Raccord: TPU

Filter/Filtre: PP

Made in China/Fabriqué en Chine

Nox Medical  
Katrínartúni 2  
IS-105 RVK  
Iceland



Caution: USA law restricts this device to sale by or on the order of a licensed medical practitioner.  
Mise en garde : La loi fédérale des États-Unis limite la vente de ce dispositif à ou à la demande d'un médecin.

NOX-FTC50-LABEL1-REV10

YYYY is replaced with the manufacturing year  
YYMMDD is replaced with the manufacturing date

## Specifications

- Storage
  - Temperature: -20°C to +50°C (0°F to 120°F)
  - Rel. Humidity: 10-95% (non-condensing)
- Shelf Life: 5 years



Nox Medical  
Katrínartúni 2 IS-105 Reykjavík Iceland  
[noxmedical@noxmedical.com](mailto:noxmedical@noxmedical.com)  
tel:+354-5707170

Technical drawing for incoming inspection

Product REF code:

552110

Part Number:

NOX-FTC50

Date:

30.07.2013

Rev:

07

Page:

1 of 2

# Shipment Packaging Instructions

- The FTC50 packaging should be packed in a larger shipment carton box.
- The box needs to be able to hold however many FTC50 are being packed in the shipment box.
- Bubble wrap should be used to fill up any empty place there might be in the shipment box.
- After the box has been closed the FTC50 should not be able to move around in it.
- One or more FTC50 can be placed in each box.
- Whenever possible, place all FTC50 being shipped in one box.



Nox Medical  
Katrínartúni 2 IS-105 Reykjavík Iceland  
noxmedical@noxmedical.com  
tel:+354-5707170

Technical drawing for incoming inspection

Product REF code:

552110

Part Number:

**NOX-FTC50**

Date:

30.07.2013

Rev:

07

Page:

2 of 2