

Nox Snap-on Electrode Cables

Nox Medical proudly announce that the
Nox Snap-On Electrode Cables have been released

The **Nox Snap-On Electrode Cables** were designed for increased ergonomics when performing sleep studies.

The electrode cables come in different lengths and colors to ease hook-up and minimize cable entanglement.

The electrode cables come in adult and pediatric lengths to suit all patients.

The electrode cables are intended to enable recording of EMG/ECG signals and to be used with industry standard stud electrode pads.

- **Unipolar Nox Snap-On Electrode Cables**

- Validated for use with the Nox T3 and Nox A1 Systems
- White colored leads are in ideal lengths for Bruxism recordings
- Green colored leads are ideal for the patient ground
- Grey and black colored leads are ideal for LM recordings
- Orange colored leads are ideal for ECG recordings
- The leads with beige colored heads are ideal for pediatrics



- **Bipolar Nox Snap-On Electrode Cables**

- Validated for use with the Nox A1 System
- The lead is divided into two separated leads
- Orange colored cables are ideal for ECG recordings
- Grey and black colored cables are ideal for LM recordings
- White colored cables are ideal for Bruxism recordings
- The leads with beige colored heads are ideal for pediatrics



Nox Snap-on Electrode Cables

Bipolar Nox Snap-on Electrode Cables, Adult and Pediatric lengths

55 43 10	50/100 cm (20/40 in), orange, keyhole connector, 1pk
55 43 11	30/50 cm (12/20 in), beige-orange, keyhole connector, 1pk
55 43 12	148/150 cm (58/60 in), grey, keyhole connector, 1pk
55 43 13	98/100 cm (39/40 in), beige-grey, keyhole connector, 1pk
55 43 14	148/150 cm (58/60 in), black, keyhole connector, 1pk
55 43 15	98/100 cm (39/40 in), beige-black, keyhole connector, 1pk
55 43 16	50/52 cm (20/21 in), white, keyhole connector, 1pk
55 43 17	30/32 cm (12/13 in), beige-white, keyhole connector, 1pk

Unipolar Nox Snap-on Electrode Cables, Adult and Pediatric lengths

55 40 20	50 cm (20 in), white, 1.5mm connector, 2pk
55 40 21	30 cm (12 in), beige-white, 1.5 mm connector, 2pk
55 40 22	100 cm (40 in), green, 1.5 mm connector, 1pk
55 40 23	50 cm (20 in), beige-green, 1.5 mm connector, 1pk
55 40 24	150 cm (60 in), grey, 1.5 mm connector, 2pk
55 40 25	100 cm (40 in), beige-grey, 1.5 mm connector, 2pk
55 40 26	150 cm (60 in), black, 1.5 mm connector, 2pk
55 40 27	100 cm (40 in), beige-black, 1.5 mm connector, 2pk
55 40 28	100 cm (40 in), orange, 1.5 mm connector, 2pk

On our [Nox Support Site](#) you can find [Technical Drawings](#) for the Nox Snap-on Electrode Cables.

Bipolar Nox Snap-on Electrode Cables, Adult and Pediatric lengths

55 43 10	NOX-ESNAP1PAC1-Technical Drawing.pdf
55 43 11	NOX-PESNAP1PAC1-Technical Drawing.pdf
55 43 12	NOX-ESNAP2PAC1-Technical Drawing.pdf
55 43 13	NOX-PESNAP2PAC1-Technical Drawing..pdf
55 43 14	NOX-ESNAP3PAC1-Technical Drawing.pdf
55 43 15	NOX-PESNAP3PAC1-Technical Drawing.pdf
55 43 16	NOX-ESNAP4PAC1-Technical Drawing.pdf
55 43 17	NOX-PESNAP4PAC1-Technical Drawing.pdf

Unipolar Nox Snap-on Electrode Cables, Adult and Pediatric lengths

55 40 20	NOX-ESNAP7PAC2-Technical Drawing.pdf
55 40 21	NOX-PESNAP7PAC2-Technical Drawing.pdf
55 40 22	NOX-ESNAP8PAC1-Technical Drawing.pdf
55 40 23	NOX-PESNAP8PAC1-Technical Drawing.pdf
55 40 24	NOX-ESNAP9PAC2-Technical Drawing.pdf
55 40 25	NOX-PESNAP9PAC2-Technical Drawing.pdf
55 40 26	NOX-ESNAP10PAC2-Technical Drawing.pdf
55 40 27	NOX-PESNAP10PAC2-Technical Drawing.pdf
55 40 28	NOX-ESNAP11PAC2-Technical Drawing.pdf

At the same time the following products will be obsolete:

Snap On Lead Wires	
55 40 12	Snap on Lead Wire, 60cm (24")
55 40 13	Snap on Lead Wire, 150cm (60")