



Nox A1 EEG Head Cable

Nox A1 EEG 5 Lead Gold Electrode Cable





Nox Medical proudly
announce that the
Nox A1 EEG Head Cable
and
**Nox A1 EEG 5 Lead Gold
Electrode Cable**
have been released

The cables were designed for
increased ergonomics when
performing sleep studies.

Nox A1 EEG Head Cable, Adult, Ref: 56 20 20

Length 90cm

Nox A1 EEG 5 Head Cable, Pediatric, Ref: 56 20 21

Length 70cm

The EEG Head Cable is anchored to the patient's forehead with a snap-on electrode. This, coupled with the unique design of the EEG Electrode Cables decreases the likelihood of the electrodes falling off during the night. This feature increases robustness of performing sleep studies.

Two Nox A1 EEG 5 Lead Gold Electrode Cables connect to the cable via mini USB.



Nox A1 EEG 5 Lead Gold Electrode Cable, Ref 54 44 11

Five different colored EEG leads with gold plated cup electrodes.

The colors of the leads are:

- Green (EOG)
- Blue (F3/F4)
- Yellow (C3/C4)
- Grey (O1/O2)
- Red (A1/A2)

The gold cup electrodes are intended to enable recording of EEG/EOG signals during sleep studies.

10mm diameter gold plated copper cup electrodes.

The electrodes are intended for both pediatric and adult patients.



On our [Nox Support Site](#) you can find [Technical Drawings](#) for incoming products inspection of the Nox A1 EEG Cables. The drawings give you information for each product and instructions for the shipment packaging.