

Product Information

NOX Abdomen Cable

- 1 unit/unité

NOX-ATC1-0141




LOT YYYMMDD
Barcode



For distributor information go to: www.noxmedical.com
Pour plus d'informations sur le distributeur : www.noxmedical.com

Made in China/Fabriqué en Chine



Caution: USA law restricts this device to sale by or on the order of a licensed medical practitioner.

Mise en garde : La loi fédérale des États-Unis limite la vente de ce dispositif à ou à la demande d'un médecin.

NOX-ATC1-LABEL1-REV8

YYYY is replaced by the actual manufacturing year

XXXXXX is replaced with the actual lot number

Barcode and barcode text is replaced with actual lot number



Specifications

- Storage
 - Temperature: -20°C to +50°C (0°F to 120°F)
 - Rel. Humidity: 10-95% (non-condensing)
- Shelf Life: 5 years



Nox Medical
Katrinartuni 2 IS-105 Reykjavik Iceland
noxmedical@noxmedical.com
tel:+354-5707170

Technical drawing for incoming inspection

Product REF code:

562010

Part Number:

NOX-ATC1

Date:

13.08.2013

Rev:


03

Page:

1 of 2

Shipment Packaging Instructions

- The ATC1 cable should be packed in a carton box.
- The box needs to be able to hold however many cables are being packed in the box.
- There should be extra space for bubble wrap.
- The bubble wrap should cover the bottom and all sides of the box and be placed on top of the products before closing the box.
- After the box has been closed the cables should not be able to move around in it
- One or more ATC1 cable can be placed in each box.
- Whenever possible, place all cables being shipped in one box.

 Nox Medical			Nox Medical Katrínartúni 2 IS-105 Reykjavík Iceland noxmedical@noxmedical.com tel:+354-5707170
Technical drawing for incoming inspection			
Product REF code:		562010	Part Number:
Date:	Rev:	Page:	NOX-ATC1
13.08.2013	03	2 of 2	