








# Bill of Material -BOM

NOX-CANAF40-0110 product shall include the following:


- Nox Nasal Pressure Cannula with filter - 40 units (NOX-CANAF)
- Product label (NOX-CANAF40-LABEL1) - see label below
- Suitable packaging material

 <p><b>NOX Nasal Pressure Cannula</b> with Filter CANAF40-0120</p> <ul style="list-style-type: none"> <li>• Adult Size - 40 units</li> <li>• Nox Pressure Connector</li> <li>• Total Length: 90cm/35in</li> </ul> <p>• Taille Adulte - 40 unités • Nox Raccord de Pression • Longueur Totale : 90cm/35in</p>	<p><b>LOT</b> Unique LOT number should fit this area</p>	<p><b>REF</b>  5 5 2 0 1 0</p>	<p>Manufactured for Nox Medical/Fabriqué pour Nox Medical Katrínartúni 2 105 Reykjavík Iceland Phone +354 5707170</p>	<p>Nox Medical Katrínartúni 2 IS-105 RVK Iceland</p> 
<p>For distributor information go to: <a href="http://www.noxmedical.com">www.noxmedical.com</a> Pour Plus d'informations sur le distributeur : <a href="http://www.noxmedical.com">www.noxmedical.com</a></p>	 YYYYY			
<p>Caution: USA law restricts this device to sale by or on the order of a licensed medical practitioner. Mise en garde : La loi fédérale des États-Unis limite la vente de ce dispositif à ou à la demande d'un médecin.</p>	<p>Assembled in Mexico/Assemblé au Mexique</p>	<p>NOX-CANAF40-LABEL1-REV5</p>		

YYYY is replaced by the actual production year.  
Insert unique LOT number in the designated area.

## Specifications

- Storage
  - Temperature: -20°C to +50°C (0°F to 120°F)
  - Rel. Humidity: 10-95% (non-condensing)
- Shelf Life: 2 years


		<p>Nox Medical Katrínartúni 2 IS-105 Reykjavík Iceland <a href="mailto:noxmedical@noxmedical.com">noxmedical@noxmedical.com</a> tel:+354-5707170</p>
<p>Technical drawing for incoming inspection</p>		
<p>Product REF code: <b>552010</b></p>		<p>Part Number:</p>
<p>Date: 14.08.2013</p>	<p>Rev: 04</p>	<p>Page: 1 of 2</p>
<p><b>NOX-CANAF40</b></p>		

# Packaging Specifications

NOX-CANAF40 product shall be packed in a suitable sized carton box to protect the product in shipping and storing environment. The box should be as small as possible to fit the 40 cannulas. The carton box should be closed with tape and packaging label should be secured to the carton box.

- The package shall be placed in a larger shipment carton for shipment.
- One or more packages can be packed together in one shipment box.
- The shipment box shall be chosen such that the number of packages being shipped fits exactly in it such that it takes as little volume as possible during shipment.



		Nox Medical Katrinartuni 2 IS-105 Reykjavik Iceland <a href="mailto:noxmedical@noxmedical.com">noxmedical@noxmedical.com</a> tel:+354-5707170
Technical drawing for incoming inspection		
Product REF code:		Part Number:
<span style="font-size: 24pt;">552010</span>		<span style="font-size: 24pt;">NOX-CANAF40</span>
Date: 14.08.2013	Rev: 04	Page: 2 of 2