

# Nox Filter Tube Connector

Ref No 55 21 10



- The **Nox Filter Tube Connector** is used to connect the Nox Cannula with Luer-lock or a mask pressure tube to the pressure lock of the Nox T3 Device.
- The package of Nox Filter Tube Connectors contains a total of 50 units. Each connector is for single use/disposable.
- The new Nox Filter Tube Connector has already been used for the Nox Cannula with Filter, ref. no. 55 20 10, and has proven to be more reliable than the former model. Therefore, the new Nox Filter Tube Connector replaces the former model.
- The package of **Nox Filter Tube Connector** is released from version 0200.

The reference number for the Filter Tube Connector package is 55 21 10.

On our [Nox Support Site](#) you can find the [Technical Drawing of Filter Tube Connector \(FTC50\)](#) and more information about Nox products.