

Product Information

NOX RIP Belts Disposable

Size Extra Large - 14 sets/Taille Extra Large - 14 paires

NOX-RIPDXL14-0111

LOT YYYMMDD
barcode

REF 
5 5 1 0 5 0

 Nox Medical

For distributor information go to: www.noxmedical.com/ Pour plus d'informations sur le distributeur : www.noxmedical.com/

Patented Design/Design breveté

Made in China/Fabriqué en Chine

 YYY



Nox Medical
Katrínartúni 2
IS-105 RVK
Iceland



Caution: USA law restricts this device to sale by or on the order of a licensed medical practitioner.


Mise en garde : La loi fédérale des États-Unis limite la vente de ce dispositif à ou à la demande d'un médecin.


NOX-RIPDXL14-LABEL1-REV1

YYYY is replaced by the actual manufacturing year
YYMMDD is replaced by the manufacturing date

Specifications

- Storage
 - Temperature: -20°C to +50°C (0°F to 120°F)
 - Rel. Humidity: 10-95% (non-condensing)
- Shelf Life: 2 years



 Nox Medical		Nox Medical Katrínartúni 2 IS-105 Reykjavík Iceland noxmedical@noxmedical.com tel: +354-5707170	
Technical drawing for incoming inspection			
Product REF code:		Part Number:	
551050		NOX-RIPDXL14	
Date:	Rev:	Page:	
15.04.2013	01	1 of 2	

Shipment Packaging Instructions



- The RIPDXL14 package should be packed in a carton box for shipping.
- The box needs to be able to hold however many packages are being packed in the box.
- The box should be selected such the RIPDXL14 packages can fill up the entire box. If there is any empty space it should be filled up with bubble wrap or carton material.
- After the box has been closed the RIPDXL14 packages should not be able to move around in it
- One or more RIPDXL14 can be placed in each box.
- RIPDXL14 packages can be packed in a box with other products from the same line such as RIPDL20, RIPDM20, RIPDS20 or RIPDP20.

 Nox Medical

Nox Medical
Katrínartúni 2 IS-105 Reykjavík Iceland
noxmedical@noxmedical.com
tel:+354-5707170

Technical drawing for incoming inspection

Product REF code:

551050

Part Number:

Date:

15.04.2013

Rev:

01

Page:

2 of 2

NOX-RIPDXL14