

New Products

Nox Gold Cup Electrode Leads

Ref No 55 44 10

Nox Medical proudly announce that the
Nox Gold Cup Electrode Leads
have been released

Nox Gold Cup Electrode Leads, Ref: 55 44 10

Length 100cm

The package includes 10 mixed color EEG leads attached with gold plated copper cup electrodes at the one end and a Touch Proof connector at the other end.

The gold cup electrodes are intended to enable recording of EEG/EOG signals during sleep studies.

Intended for both pediatric and adult patients.



On our [Nox Support Site](#) you can find a [Technical Drawing](#) for incoming products inspection of the Nox Gold Cup Electrode Leads. The drawing gives you information for the product and instructions for the shipment packaging.