

Product Specifications



Bill of Material -BOM

NOX-CANA50-0100 product shall include the following:

- Nox Nasal Pressure Cannula with Luer-Lock - 50 units (NOX-CANA-0100)
- Product label (NOX-CANA50-LABEL1) - see label below
- Brown carton box 23,5x23,5x16,5cm

 NOX Nasal Pressure Cannula with Luer-Lock CANA50-0100	LOT Unique LOT number should fit this area
<ul style="list-style-type: none">• Adult Size - 50 units• Male Luer Connector• Total Length: 90cm/35in	REF 5 5 2 0 2 0
<p>For distributor information go to: www.noxmedical.com Pour Plus d'informations sur le distributeur : www.noxmedical.com</p> <p>Caution: USA law restricts this device to sale by or on the order of a licensed medical practitioner. Mise en garde : La loi fédérale des États-Unis limite la vente de ce dispositif à ou à la demande d'un médecin.</p>	<p>Manufactured for Nox Medical/Fabriqué pour Nox Medical Katrínartúni 2 105 Reykjavík Iceland Phone +354 5707170</p> <p>Nox Medical Katrínartúni 2 IS-105 RVK Iceland</p> <p> -20 °C (0°F) to +50 °C (120°F)</p> <p> YYYYY</p> <p>Assembled in Mexico/Assemblé au Mexique</p> <p>NOX-CANA50-LABEL1-REV4</p>

Comments:

Replace YYYYY with the manufacturing year.
Insert LOT number in the designated area.

		Nox Medical Katrínartúni 2 IS-105 Reykjavík Iceland noxmedical@noxmedical.com tel:+354-5707170
Technical drawing for incoming inspection		
Product REF code: 552020		Part Number:
Date: 13.08.2013	Rev: 05	Page: 1 of 2
		NOX-CANA50

Packaging Specifications

NOX-CANA50-0100 product shall be packed in a 23,5x23,5x16,5 cm carton box to protect the product in shipping and storing environment.

The carton box should be closed with tape and packaging label should be secured to the carton box.



Nox Medical

Nox Medical
 Katrinartuni 2 IS-105 Reykjavik Iceland
norxmedical@noxmedical.com
 tel:+354-5707170

Technical drawing for incoming inspection

Product REF code:		552020		Part Number:	
Date:	Rev:	Page:	NOX-CANA50		
13.08.2013	05	2 of 2			