

nox medical

NOX A1 PSG SYSTEM

With Noxturnal



Hardware	Comment	Item codes	PSG ambulatory	PSG online basic (w/o video)	PSG online regular w/ video	PSG online extended with all DC and w/video
- Nox A1 System with SPO2, CE	- Complete HW Nox A1 System with SPO2 (sufficient for ambulatory use)	51 30 10	x	x	x	x
- C1 Lab Control Unit	- Required HW to receive A1 online streaming to LAN	54 40 20		x	x	x
Software and Tablet Application						
- Noxturnal A1 PSG File Format License	- Nox A1 preparation / download ambulatory A1 recordings.	53 40 10	x			
- Noxturnal A1 Online Acquisition License	- Nox A1 online acquisition without video, as well as Nox A1 preparation / download ambulatory A1 recordings, reading Nox A1 format recordings incl. video playback	53 60 10		x	x	x
- Noxturnal Video Online Acquisition License	- This enables online video acquisition	53 60 11			x	x
- Noxturnal DC Channel Online Acquisition License	- This enables online DC channel acquisition from C1 (applied per C1 device)	53 60 13				x
Sensors and System Extensions						
- Nox A1 Starter Sensor Kit, Adult		55 84 10	x	x	x	x
Optional items / functionality						
- Nox A1 Starter Sensor Kit, Pediatric		55 84 11	Option	Option	Option	Option
- Nox W7 Link - S	- Bluetooth connection for Sentec - Capnography in ambulatory use (Sentec: SDM)	54 40 10	Option	Not required	Not required	Not required
- Nox W7 Link - R	- Bluetooth connection for AUX devices type R in ambulatory use (Radiometer: Tosca TCM 40, CombiM TCM4; ResMed: ResMed Airsense)	54 40 11	Option	Not required	Not required	Not required
- Nox W7 Link - A	- Bluetooth connection for AUX devices type Airsense in ambulatory use	54 40 12	Option	Not required	Not required	Not required
- Android tablet	- Regular Android tablets / phones should be compatible (few exception may apply).		Rec.	Option	Option	Option
- Video camera for online use	- Nox offers specifications and model recommendations				Required	Required
- Network hub for additional flexibility in online use	- Network hub to gather multiple IP devices in online studies within the patient room.			Option	Option	Option



	PSG ambulatory	PSG online basic	PSG online basic w/video	PSG online extended
Nox offers	- NOX A1 PSG System - SW PSG A1 (ambulat.) - SW for tablet - SSK (40 studies)	- NOX A1 PSG System - HW for online link - SW PSG online - SW for tablet - SSK (40 studies)	- NOX A1 PSG System - HW for online link - SW online and video - SW for tablet - SSK (40 studies)	- NOX A1 PSG System - Nox C1 Access Point - SW online and video - SW for tablet - SSK (40 studies)
Distributor adds	- Tablet	- Tablet - Network Hub	- Tablet - Network Hub - (IP) Video System	- Tablet - Network Hub - (IP) Video System
Existing resources	- Computer and monitor - Soundbar (snoring) - Cables	- Computer and monitor - Soundbar (snoring) - Intercom - Cables	- Computer and monitor - Soundbar (snoring) - Intercom - Cables	- Computer and monitor - Soundbar (snoring) - Intercom - Cables

